

## Deltalab 2ml Polypropylene Graduated Microtubes

**Product Code:** 4092.7N

**Brand:** Deltalab

Made of ultra clear polypropylene, Deltalab's 2ml Polypropylene Graduated Microtubes have a frosted area for writing.

Their flat cap has a secure lock with a safety shield.

These microtubes can be handled easily with just one hand.

They can be used at temperatures from -20°C up to 121°C and can withstand centrifugation up to 12,000xg.

### Specifications

<b>Capacity (ml)</b>	2
<b>Dimensions (Ø x h)</b>	10.3 x 40.5mm
<b>Quantity</b>	5000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)