

Deltalab Graduated Beakers (250ml)

Product Code: 409250

Brand: Deltalab

Ideal for general laboratory use, Deltalab's Graduated Beakers (250ml) are made of durable polypropylene for safe use. These lab beakers are autoclavable up to 121°C, have clearly printed graduations, and are in compliance with ISO 6706 and BS 5404 standards.



Specifications

| | |
|----------------------|---------------|
| Capacity (ml) | 250 |
| Diameter (mm) | 71 |
| Height (mm) | 96 |
| Material | Polypropylene |
| Sterility | Non-Sterile |
| Quantity | 16 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com