

Medline Polypropylene Clear Dipper Samplers with Handles (30ml)

Product Code: 409277

Brand: Medline Scientific

Medline's Polypropylene Clear Dipper Samplers with Handles (30ml) allow laboratories to carry out simple sample collections and subsequent transportation. The handle neatly snaps off after sample collection from the vat for the ease of transportation.

Each dipper is supplied with a flow-seal screw cap to ensure leak-proof closure.

Specifications

Colour	Clear
Volume (ml)	30
Handle Length (mm)	385
Material	Polypropylene
Inner Pack	1
Quantity	50

