

Medline Polypropylene Blue Dipper Samplers with Handles (30ml)

Product Code: 409278

Brand: Medline Scientific

Medline's Polypropylene Blue Dipper Samplers with Handles (30ml) allow laboratories to carry out simple sample collections and subsequent transportation. The handle neatly snaps off after sample collection from the vat for the ease of transportation.

Each dipper is supplied with a flow-seal screw cap to ensure leak-proof closure.

Specifications

Colour	Blue
Volume (ml)	30
Handle Length (mm)	385
Material	Polypropylene
Inner Pack	1
Quantity	50

