


## Deltalab 0.2ml Blue PCR Tubes

**Product Code:** 4094.1A

**Brand:** Deltalab

Made of polypropylene, Deltalab's 0.2ml Blue PCR Tubes feature attached hinged caps. The caps are flat and easily pierceable.

Autoclavable, these PCR tubes are certified DNase, RNase, and PCR inhibitors free.



### Specifications

<b>Colour</b>	Blue
<b>Quantity</b>	1000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)