

Deltalab 0.2ml Yellow PCR Tubes

Product Code: 4094.1A

Brand: Deltalab

Made of polypropylene, Deltalab's 0.2ml Yellow PCR Tubes feature attached hinged caps. The caps are flat and easily pierceable.

Autoclavable, these PCR tubes are certified DNase, RNase, and PCR inhibitors free.



Specifications

Colour	Yellow
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com