

Deltalab 0.2ml Red PCR Tubes

Product Code: 4094.1R


Brand: Deltalab

Made of polypropylene, Deltalab's 0.2ml Natural PCR Tubes feature attached hinged caps. The caps are flat and easily pierceable.

Autoclavable, these PCR tubes are certified DNase, RNase, and PCR inhibitors free.

Specifications

Colour	Red
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com