

## Deltalab 0.5ml Blue Graduated PCR Tubes

**Product Code:** 4094.2A


**Brand:** Deltalab

Made of autoclavable polypropylene, Deltalab's 0.5ml Blue Graduated PCR Tubes have attached hinged caps, which are flat and easily pierceable.

The tubes can be easily opened and closed with one hand.

They have moulded-in graduations in 0.1 ml increments from 0.1 up to 0.6ml, and a fronted panel on their side for writing or labelling.

These PCR tubes are certified DNase, RNase, and PCR inhibitors free.



### Specifications

<b>Colour</b>	Blue
<b>Quantity</b>	1000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)