

Deltalab 0.5ml Natural Graduated PCR Tubes

Product Code: 4094.2N

Brand: Deltalab

Made of autoclavable polypropylene, Deltalab's 0.5ml Natural Graduated PCR Tubes have attached hinged caps, which are flat and easily pierceable.

The tubes can be easily opened and closed with one hand.

They have moulded-in graduations in 0.1 ml increments from 0.1 up to 0.6ml, and a fronted panel on their side for writing or labelling.

These PCR tubes are certified DNase, RNase, and PCR inhibitors free.



Specifications

Colour	Natural
Quantity	1000

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com