

Deltalab 0.2ml PCR Tubes in Strips (Low Profile)

Product Code: 4094.4N

Brand: Deltalab

Made of polypropylene, Deltalab's 0.2ml PCR Tubes in Strips (Low Profile) have attached hinged caps. They are autoclavable and suitable for real time PCR.


The caps are flat and easily pierceable.

They have moulded-in graduations in 0.1 ml increments from 0.1 up to 0.6ml, and a fronted panel on their side for writing or labelling.

These PCR tubes are certified DNase, RNase, and PCR inhibitors free.

Specifications

Description	Strip of 8 tubes and 8 attached caps, low profile
Quantity	120



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com