

## Deltalab 0.2ml Real Time PCR Tubes in Strips

**Product Code:** 4094.5BP

**Brand:** Deltalab

Made of polypropylene, Deltalab's 0.2ml Real Time PCR Tubes in Strips are autoclavable and feature attached flat caps.

The strips are manufactured by biomaterial moulding so the tubes are made of opaque white polypropylene, while the caps are compounded of transparent, optically clear polypropylene.


These opaque white tubes perform the highest performance in real time PCR as they avoid crosstalk between wells.

The caps are easily pierceable.

These PCR tubes are certified DNase, RNase, and PCR inhibitors free. They come in bags of 10 strips.

### Specifications

<b>Description</b>	Strip of 8 white 0.2ml QPCR tubes and 8 attached caps
<b>Quantity</b>	120



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)