

Deltalab 0.2ml Real Time PCR Tubes

Product Code: 4094.5N

Brand: Deltalab

Made of polypropylene, Deltalab's 0.2ml Real Time PCR Tubes are autoclavable and feature attached hinged caps.

The caps are easily pierceable and offer optical quality, thus allowing their application in real time PCR.

These PCR tubes are certified DNase, RNase, and PCR inhibitors free.

Specifications

Description	PCR 0.2ml QPCR tube
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com