

Deltalab 40ml Polypropylene Containers (30 x 70mm), Individually Wrapped

Product Code: 409426

Brand: Deltalab

With polyethylene caps, Deltalab's Polypropylene 40ml Containers (30 x 70mm), Individually Wrapped are ideal for liquid specimen collection, transportation, and storage in rough conditions. These graduated sample containers have a strong gripping, which means they can be opened and closed easily when wearing gloves. Their clear graduations enable simple viewing.

Specifications

Volume (ml)	40
Max. Volume (ml)	39
Recommended Volume (ml)	30
Description	Ultra clear polypropylene, individually wrapped
Dimensions (mm)	30 x 70
Quantity	800



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com