


Deltalab 8-Strip 0.2ml PCR Tubes (12 Tube Strip Attached to 12 Domed Cap Strips)

Product Code: 4095.4N

Brand: Deltalab

Made of natural colour polypropylene, Deltalab's 0.2ml PCR Tubes (12 Tube Strip Attached to 12 Domed Cap Strips) are autoclavable and certified DNase, RNase, and PCR inhibitors free.

The cap and tube strips allow easy handling and inventory of tubes and caps. The cap strips prevent cross-contamination from tube to tube.



Specifications

Description	12 tube strip attached to 12 domed cap strips
Quantity	80

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com