


Deltalab 8-Strip 0.2ml PCR Tubes (Strip of 8 Flat Caps)

Product Code: 4095.7N

Brand: Deltalab

Made of natural colour polypropylene, Deltalab's 0.2ml PCR Tubes (Strip of 8 Flat Caps) are autoclavable and certified DNase, RNase, and PCR inhibitors free.

These cap strips allow easy handling and inventory of caps. They also prevent cross-contamination from tube to tube.



Specifications

Description	Strip of 8 flat caps
Quantity	125

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com