

Deltalab 60ml Aseptic Polypropylene Containers (38 x 65mm), Labelled

Product Code: 409501

Brand: Deltalab

Part of a new generation of sterile sample containers with high safety closure biflex caps, Deltalab's 60ml Aseptic Polypropylene Containers (38 x 65mm), Labelled are specially designed for pneumatic transport of specimens.



Specifications

| | |
|--------------------------------|--------------------------|
| Volume (ml) | 60 |
| Max. Volume (ml) | 62 |
| Recommended Volume (ml) | 50 |
| Description | Polypropylene (labelled) |
| Dimensions (mm) | 38 x 65 |
| Quantity | 600 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com