

Deltalab 60ml Polypropylene Containers (38 x 65mm), Individually Wrapped

Product Code: 409526G

Brand: Deltalab

Deltalab's 60ml Polypropylene Containers (38 x 65mm), Individually Wrapped have moulded-in graduations in 5ml increments.

Each laboratory container comes with a red polyethylene biflex cap. They are supplied capped.

They are ideal for pneumatic transport because they have a watertight closing system.

Specifications

Description	Polypropylene, individually wrapped
Sterility	Aseptic
Quantity	600



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com