


## Deltalab 60ml Aseptic Polypropylene Containers (38 x 65mm) With Spoon

**Product Code:** 409552

**Brand:** Deltalab

Part of a new generation of sterile sample containers with high safety closure biflex caps, Deltalab's 60ml Aseptic Polypropylene Containers (38 x 65mm) With Spoon are specially designed for pneumatic transport of specimens.



### Specifications

<b>Volume (ml)</b>	60
<b>Max. Volume (ml)</b>	62
<b>Recommended Volume (ml)</b>	50
<b>Description</b>	Polypropylene with Spoon
<b>Dimensions (mm)</b>	38 x 65
<b>Quantity</b>	600

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)