

## Deltalab 0.2ml Natural PCR Tubes in Strips (With Flat Caps)

**Product Code:** 4096.2N

**Brand:** Deltalab

Made of autoclavable polypropylene, Deltalab's 0.2ml Natural PCR Tubes in Strips (With Flat Caps) are certified RNase, DNase, and pyrogen free.

The strips include 8 integral 0.2ml tubes with ultra thin sidewalls and bottoms with individually attached caps. Every cap embodies a shield in order to prevent contamination when being opened. While easily opened and sealed with one hand, their positive sealing will fully protect contents from evaporation during the whole thermal cycle. These strips can be cut to any length while each tube has its own attached cap.

### Specifications

Cap	Flat cap
Colour	Natural
Quantity	125



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)