


Deltalab 0.2ml Natural PCR Tubes in Strips (With Domed Caps)

Product Code: 4096.3N

Brand: Deltalab

Made of autoclavable polypropylene, Deltalab's 0.2ml Natural PCR Tubes in Strips (With Domed Caps) are certified RNase, DNase, and pyrogen free.

The strips include 8 integral 0.2ml tubes with ultra thin sidewalls and bottoms with individually attached caps. Every cap embodies a shield in order to prevent contamination when being opened. While easily opened and sealed with one hand, their positive sealing will fully protect contents from evaporation during the whole thermal cycle. These strips can be cut to any length while each tube has its own attached cap.



Specifications

Cap	Domed cap
Colour	Natural
Quantity	125

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com