


Deltalab 30ml Non-Sterile Polypropylene Universal Containers (25 x 90mm) With Polystyrene Spoon

Product Code: 409602

Brand: Deltalab

Supplied non-sterile and autoclavable, Deltalab's Non-Sterile 30ml Polypropylene Universal Containers (25 x 90mm) With PS Spoon are manufactured with an innovative process resulting in one of the most resistant containers available on the market. The screw closure can assure an almost leak-proof closing system.



Specifications

Volume (ml)	30
Max. Volume (ml)	31.1
Recommended Volume (ml)	25
Description	Polypropylene with PS spoon
Dimensions (mm)	25 x 90
Sterility	Non-Sterile
Quantity	400

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com