

## Deltalab 150ml Aseptic Polypropylene Containers (57 x 73mm), Labelled

**Product Code:** 409701

**Brand:** Deltalab

Manufactured with ultra clear polypropylene, Deltalab's 150ml Aseptic Polypropylene Containers (57 x 73mm), Labelled feature a polyethylene leak proof red screw cap and a frosted panel on the side for easy writing.

Each sample container has moulded graduations in 20ml increments up to 100ml.

### Specifications

<b>Max. Volume (ml)</b>	150
<b>Recommended Volume (ml)</b>	125
<b>Description</b>	Polypropylene labelled
<b>Sterility</b>	Aseptic
<b>Quantity</b>	450

