

Deltalab 150ml Uncapped Polypropylene Containers (57 x 73mm)

Product Code: 409703

Brand: Deltalab

Manufactured with ultra clear polypropylene, Deltalab's 150ml Uncapped Polypropylene Containers (57 x 73mm) feature polyethylene leak proof red screw caps supplied in the same box as the containers and a frosted panel on the side for easy writing.

Each sample container has moulded graduations in 20ml increments up to 100ml.

Specifications

Max. Volume (ml)	150
Recommended Volume (ml)	125
Description	Polypropylene uncapped
Quantity	620

