

## Dotalab 150ml Sterile R Polypropylene Containers (57 x 73mm), Individually Wrapped

**Product Code:** 409726.G

**Brand:** Dotalab

Manufactured with ultra clear polypropylene, Dotalab's 150ml Sterile R Polypropylene Containers (57 x 73mm), Individually Wrapped feature a polyethylene leak proof red screw cap and a frosted panel on the side for easy writing.

Each sample container has moulded graduations in 20ml increments up to 100ml.

### Specifications

<b>Max. Volume (ml)</b>	150
<b>Recommended Volume (ml)</b>	125
<b>Description</b>	Polypropylene individually wrapped
<b>Sterility</b>	Sterile R
<b>Quantity</b>	350

