

Deltalab 150ml Sterile A Polypropylene Containers (57 x 73mm), Individually Wrapped

Product Code: 409726

Brand: Deltalab

Manufactured with ultra clear polypropylene, Deltalab's 150ml Sterile A Polypropylene Containers (57 x 73mm), Individually Wrapped feature a polyethylene leak proof red screw cap and a frosted panel on the side for easy writing.

Each sample container has moulded graduations in 20ml increments up to 100ml.

Specifications

Max. Volume (ml)	150
Recommended Volume (ml)	125
Description	Polypropylene individually wrapped
Sterility	Sterile A
Quantity	350

