

## Deltalab 150ml Non-Sterile Polypropylene Containers with Spoon (57 x 73mm)

**Product Code:** 409752

**Brand:** Deltalab

Ideal for sample collection and transport, Deltalab's 150ml Non-Sterile Polypropylene Containers with Spoon (57 x 73mm) are manufactured with polyethylene leak proof screw caps and polystyrene spoons.

Graduated with 100ml increments, these sample containers also feature a frosted panel on the side for easy writing.



### Specifications

<b>Max. Volume (ml)</b>	150
<b>Recommended Volume (ml)</b>	125
<b>Description</b>	Polypropylene with PS spoon
<b>Dimensions (mm)</b>	57 x 73
<b>Sterility</b>	Non-Sterile
<b>Quantity</b>	450

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)