

Deltalab 150ml Non-Sterile Polypropylene Containers with Spoon (57 x 73mm), Individually Wrapped

Product Code: 409756

Brand: Deltalab

Ideal for sample collection and transport, Deltalab's 150ml Non-Sterile Polypropylene Containers with Spoon (57 x 73mm), Individually Wrapped are manufactured with polyethylene leak proof screw caps and polystyrene spoons.

Graduated with 100ml increments, these sample containers also feature a frosted panel on the side for easy writing.



Specifications

Max. Volume (ml)	150
Recommended Volume (ml)	125
Description	Polypropylene with PS spoon, individually wrapped
Dimensions (mm)	57 x 73
Sterility	Non-Sterile
Quantity	350

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com