

## Sterile Conical High Resistance Tubes (15ml)

**Product Code:** 409920

**Brand:** Medline

Made of transparent and copolymer polypropylene, these 15ml Sterile Conical Tubes are autoclavable and feature a solvent resistant white panel and black graduations for both dark or clear samples.

Their blue caps are made of polyethylene with an internal liner (the liner is not autoclavable), which ensures leakproof integrity.

Each tube is sterilised by radiation and withstands temperatures down to -90°C.



### Specifications

<b>Capacity (ml)</b>	15
<b>Dimensions (mm)</b>	17 x 118
<b>Resistance to Centrifugation</b>	17,000xg
<b>Material(s)</b>	Polypropylene (tube) and Polyethylene (cap)
<b>Sterility</b>	Sterile
<b>Quantity</b>	500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)