

Medline 50ml Racked Conical High-Resistance

Tubes Product Code: 409922

Brand: Medline Scientific

Medline's 50ml Conical High-Resistance Tubes are made of transparent, copolymer polypropylene. The caps are made of polyethylene with an internal liner, which ensures a leakproof seal.

They are DNase, RNase, endotoxins, and metal free.

They are sterile by radiation and can withstand temperatures down to -90°C (15ml) and -80°C (50ml).

These tubes are also autoclavable (the liner of the caps are not autoclavable).



Specifications

Capacity (ml)	50
Dimensions (mm)	29.6 x 114.6
Material(s)	Polypropylene (tubes) and Polyethylene (caps)
Sterility	Sterile
Quantity	10 x 50

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com