

## Deltalab 50ml Hinged Lid Containers

**Product Code:** 410046

**Brand:** Deltalab

Deltalab's 50ml Hinged Lid Containers are specially designed for milk, wine and agrolimentary oil. They can also be used for other solid or liquid samples.

These sample containers are recommended for their use with FOSS® and DELTA® analytical systems, among others. They are made of autoclavable polypropylene and feature a tamper evident lid attached to the body by means of a hinge.

The lid is conceived to maintain 90° after its opening, enabling the containers to be used with just one hand. They have a frosted square surface for identification, and have excellent mechanical and chemical resistance.



### Specifications

<b>Description</b>	50ml container, uncapped
<b>Material</b>	Polypropylene
<b>Quantity</b>	800

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)