

Deltalab 50ml Polypropylene Tamper Evident Containers (38 x 70mm), Individually Wrapped

Product Code: 419826

Brand: Deltalab

Deltalab's 50ml Polypropylene Tamper Evident Containers (38 x 70mm), Individually Wrapped feature a specially designed tamper evident ring.

These sample containers are supplied capped but not fully locked so that the container is not fully sealed.

Once samples (e.g. urine, milk, water) are collected, screw the cap in a normal way and then screw again to fully position the tamper evident ring.



Specifications

Volume (ml)	50
Max. Volume (ml)	56
Recommended Volume (ml)	45
Description	Ultra clear polypropylene, individually wrapped
Dimensions (mm)	38 x 70
Quantity	600

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com