

## Deltalab EUROTUBO 50ml Non-Skirted Conical Tubes, Unscrewed Cap

**Product Code:** 429900SP

**Brand:** Deltalab

Made of transparent polypropylene, Deltalab's EUROTUBO 50ml Non-Skirted Conical Tubes, Unscrewed Cap are suitable for centrifugation tests in various laboratories (e.g. immunology, microbiology).

Each conical tube has a continuous thread, moulded external graduations (each 5ml), and blue polyethylene caps with internal elastic obturation to assure an hermetical closure.

These tubes are autoclavable at 121°C with the cap positioned on the thread but uncapped.



### Specifications

Description	Non-Skirted, unscrewed cap
<b>Sterile</b>	No
<b>External Cap (Ø mm)</b>	34.4
<b>Internal Tube (Ø mm)</b>	27.2
<b>External Tube (Ø mm)</b>	29.5
<b>Length with Cap (mm)</b>	117.5
<b>Quantity</b>	500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)