

## Deltalab EUROTUBO 50ml Skirted Conical Tubes

**Product Code:** 429901

**Brand:** Deltalab

Made of transparent polypropylene, Deltalab's EUROTUBO 50ml Skirted Conical Tubes are suitable for centrifugation tests in various laboratories (e.g. immunology, microbiology).

Each conical tube has a continuous thread, moulded external graduations (each 5ml), and blue polyethylene caps with internal elastic obturation to assure an hermetical closure.

These tubes are autoclavable at 121°C with the cap positioned on the thread but uncapped.



### Specifications

Description	Skirted
Sterile	No
External Cap (Ø mm)	34.4
Internal Tube (Ø mm)	27.2
External Tube (Ø mm)	29.5
Length with Cap (mm)	117.5
Quantity	500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

Tel: (+44) 01865 400321

Email: [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

Website: [www.medlinescientific.com](http://www.medlinescientific.com)