

## Deltalab EUROTUBO 15ml Non-Sterile Conical Tubes

**Product Code:** 429910

**Brand:** Deltalab

Made of transparent polypropylene, Deltalab's EUROTUBO 15ml Non-Sterile Conical Tubes are suitable for centrifugation tests in various laboratories (e.g. immunology, microbiology).

Each conical tube has a continuous thread, moulded external graduations (each 0.5ml), and blue polyethylene caps with internal elastic obturation to assure an hermetical closure.

These tubes are autoclavable at 121°C with the cap positioned on the thread but uncapped.



### Specifications

| Description                 | Non-Sterile |
|-----------------------------|-------------|
| <b>Sterile</b>              | No          |
| <b>External Cap (Ø mm)</b>  | 20.9        |
| <b>Internal Tube (Ø mm)</b> | 14.35       |
| <b>External Tube (Ø mm)</b> | 17          |
| <b>Length with Cap (mm)</b> | 120         |
| <b>Quantity</b>             | 500         |

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)