

## Deltalab EUROTUBO 15ml Sterile Conical Tubes, 100 Units Per Bag

**Product Code:** 429920

**Brand:** Deltalab

Made of transparent polypropylene, Deltalab's EUROTUBO 15ml Sterile Conical Tubes, 100 Units Per Bag are suitable for centrifugation tests in various laboratories (e.g. immunology, microbiology).

Each conical tube has a continuous thread, moulded external graduations (each 0.5ml), and blue polyethylene caps with internal elastic obturation to assure an hermetical closure.

These tubes are not suitable for autoclaving.



### Specifications

<b>Description</b>	100 units in bag
<b>Sterile</b>	Sterile R
<b>External Cap (Ø mm)</b>	20.9
<b>Internal Tube (Ø mm)</b>	14.35
<b>External Tube (Ø mm)</b>	17
<b>Length with Cap (mm)</b>	120
<b>Quantity</b>	5 x 100

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)