

Deltalab EUROTUBO 15ml Sterile Conical Tubes, 100 Units Per Bag

Product Code: 429920

Brand: Deltalab

Made of transparent polypropylene, Deltalab's EUROTUBO 15ml Sterile Conical Tubes, 100 Units Per Bag are suitable for centrifugation tests in various laboratories (e.g. immunology, microbiology).

Each conical tube has a continuous thread, moulded external graduations (each 0.5ml), and blue polyethylene caps with internal elastic obturation to assure an hermetical closure.

These tubes are autoclavable at 121°C with the cap positioned on the thread but uncapped.

Specifications

Description	100 units in bag
Sterile	Sterile R
External Cap (Ø mm)	20.9
Internal Tube (Ø mm)	14.35
External Tube (Ø mm)	17
Length with Cap (mm)	120
Quantity	5 x 100

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com