

Deltalab EUROTUBO 50ml Non-Skirted Conical Tubes, Individually Wrapped

Product Code: 429926

Brand: Deltalab

Made of transparent polypropylene, Deltalab's EUROTUBO 50ml Non-Skirted Conical Tubes, Individually Wrapped are suitable for centrifugation tests in various laboratories (e.g. immunology, microbiology).

Each conical tube has a continuous thread, moulded external graduations (each 5ml), and blue polyethylene caps with internal elastic obturation to assure an hermetical closure.

These tubes are not suitable for autoclaving.



Specifications

Description	Non-Skirted, individually wrapped
Sterile	Sterile R
External Cap (Ø mm)	34.4
Internal Tube (Ø mm)	27.2
External Tube (Ø mm)	29.5
Length with Cap (mm)	117.5
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com