

Deltalab EUROTUBO 50ml Skirted Conical Tubes, Individually Wrapped

Product Code: 429927

Brand: Deltalab

Made of transparent polypropylene, Deltalab's EUROTUBO 50ml Skirted Conical Tubes, Individually Wrapped are suitable for centrifugation tests in various laboratories (e.g. immunology, microbiology).

Each conical tube has a continuous thread, moulded external graduations (each 5ml), and blue polyethylene caps with internal elastic obturation to assure an hermetical closure.

These tubes are not suitable for autoclaving.



Specifications

| Description | Skirted, individually wrapped |
|-----------------------------|-------------------------------|
| Sterile | Sterile R |
| External Cap (Ø mm) | 34.4 |
| Internal Tube (Ø mm) | 27.2 |
| External Tube (Ø mm) | 29.5 |
| Length with Cap (mm) | 117.5 |
| Quantity | 500 |

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com