

Deltalab Polypropylene Histology Moulds (24 x 24 x 5mm)

Product Code: 440132


Brand: Deltalab

The smooth interior of Deltalab's Polypropylene Histology Moulds (15 x 15 x 5mm) ensures that the paraffin does not adhere to the mould. The special material used provides excellent thermal exchange.

They are suitable for most cassettes.

Specifications

External Dimensions (mm)	50 x 37.2 x 12
Well Dimensions (mm)	24 x 24 x 5
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com