

## Deltalab White Embedding Cassettes Without Lids

**Product Code:** 440141

**Brand:** Deltalab

Made of polyacetylene, Deltalab's White Embedding Cassettes Without Lids have 34 round holes with a diameter of 2mm. They are resistant to solvents used in histological processes. It is possible to write on them.

### Specifications

<b>External Dimensions (mm)</b>	40 x 28 x 6
<b>Colour</b>	White
<b>Quantity</b>	3000



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)