

Deltalab Yellow Embedding Cassettes Without Lids


Product Code: 440143

Brand: Deltalab

Made of polyacetylene, Deltalab's Yellow Embedding Cassettes Without Lids have 34 round holes with a diameter of 2mm. They are resistant to solvents used in histological processes. It is possible to write on them.

Specifications

External Dimensions (mm)	40 x 28 x 6
Colour	Yellow
Quantity	3000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com