

## Deltalab Non-Sterile Cytology Brushes

**Product Code:** 440150

**Brand:** Deltalab

Used for cervical specimen collection, Deltalab's Non-Sterile Cytology Brushes have a blue ABS shaft and are designed for taking cell samples without damaging them.

They also include a soft tip to provide more comfort for patients.

### Specifications

<b>Total Length (cm)</b>	19.5
<b>Description</b>	Non-Sterile
<b>Quantity</b>	2500



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)