

Deltalab White Mega Cassettes

Product Code: 44125

Brand: Deltalab

Made from polyacetal, Deltalab's White Mega Cassettes have an integral lid and are designed for larger tissue samples.

Specifications

Holes (mm)	1 x 5
Cassette Dimensions (mm)	26 x 40 x 13
Colour	White
Quantity	500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com