


## Deltalab Yellow Tissue Embedding Cassettes (Histosette II) - Lids Pre-Mounted

**Product Code:** 44143

**Brand:** Deltalab

Made of acetal polymer, Deltalab's Yellow Tissue Embedding Cassettes (Histosette II) - Lid Pre-Mounted have 64 cuts in the base and lid.

The lids come pre-mounted.



### Features

- Hinged lid can be opened and closed as often as necessary
- Supplied assembled, saving time and making storage easier
- Air vents allowing more efficient filling with paraffin
- Easy handling with one hand
- Safety closure
- Writing area slanted at a 45° angle, suitable to be used with most labelling instruments

### Specifications

<b>Colour</b>	Yellow
<b>External Dimensions (mm)</b>	41 x 28 x 6
<b>Quantity</b>	1500

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)