

Deltalab Green Tissue Embedding Cassettes (Histosette II) - Lid Pre-Mounted

Product Code: 44147

Brand: Deltalab

Made of acetal polymer, Deltalab's Green Tissue Embedding Cassettes (Histosette II) - Lid Pre-Mounted have 64 cuts in the base and lid.

The lids come pre-mounted.

Features

- Hinged lid can be opened and closed as often as necessary
- Supplied assembled, saving time and making storage easier
- Air vents allowing more efficient filling with paraffin
- Easy handling with one hand
- Safety closure
- Writing area slanted at a 45° angle, suitable to be used with most labelling instruments

Specifications

Colour	Green
External Dimensions (mm)	41 x 28 x 6
Quantity	1500

Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com

