

## Deltalab Pink Microsample Cassettes with Lids

**Product Code:** 44240

**Brand:** Deltalab

Designed for small biopsies, Deltalab's Pink Microsample Cassettes with Lids have six individual cells to hold six samples. Each cell is 7 x 12mm and the lids allow multiple opening and closing.

A frosted area allows writing on the cassette.

Each cassette has near 3000 micro holes of 0.26mm diameter.

### Specifications:

<b>Colour</b>	Pink
<b>Material</b>	Polyacetylene
<b>Quantity</b>	1000

Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)