

Deltalab Blue Microsample Cassettes with Lids

Product Code: 44245

Brand: Deltalab


Designed for small biopsies, Deltalab's Blue Microsample Cassettes with Lids have six individual cells to hold six samples. Each cell is 7 x 12mm and the lids allow multiple opening and closing.

A frosted area allows writing on the cassette.

Each cassette has near 3000 micro holes of 0.26mm diameter.

Specifications:

Colour	Blue
Material	Polyacetylene
Quantity	1000



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com