

Deltalab 250ml Narrow Neck High Resistance Bottles

Product Code: 44260000

Brand: Deltalab

Designed for the storage of acids and bases, Deltalab's 250ml Narrow Neck High Resistance Bottles are made of translucent, white, high density polyethylene.

Each bottle has a blue screw cap on a combination of high and low density polyethylene. The screw caps are ribbed for a better handling.

These heavy-duty containers also have sloping shoulders for easy pouring and are supplied uncapped.

Specifications

Volume (ml)	250
Neck (Ø mm)	30
Bottom (Ø mm)	67
Height (mm)	157
Quantity	90



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com