

Deltalab 500ml Wash Bottles

Product Code: 44300000

Brand: Deltalab

Made of a combination of high and low density polyethylene, Deltalab's 250ml Wash Bottles are translucent with HDPE caps (blue), a delivery tube, and seal caps (natural).

The caps are ribbed for a better handling.

These laboratory containers are highly resistant towards all liquids, acids, and bases.

They also feature a large neck for safe and easy filling.

Specifications

Volume (ml)	500
Neck (mm)	30
Base (mm)	75
Height (mm)	185
Quantity	60



Unit 3, Tower Business Park
Warpsgrove Lane
Chalgrove
Oxfordshire
OX44 7XZ

Tel: (+44) 01865 400321

Email: enquiries@medlinescientific.com

Website: www.medlinescientific.com