

## Deltalab 1000ml Wash Bottles

**Product Code:** 44310000

**Brand:** Deltalab

Made of a combination of high and low density polyethylene, Deltalab's 1000ml Wash Bottles are translucent with HDPE caps (blue), a delivery tube, and seal caps (natural).

The caps are ribbed for a better handling.

These laboratory containers are highly resistant towards all liquids, acids, and bases.

They also feature a large neck for safe and easy filling.

### Specifications

<b>Volume (ml)</b>	1000
<b>Neck (mm)</b>	30
<b>Base (mm)</b>	90
<b>Height (mm)</b>	222
<b>Quantity</b>	45



Unit 3, Tower Business Park  
Warpsgrove Lane  
Chalgrove  
Oxfordshire  
OX44 7XZ

**Tel:** (+44) 01865 400321

**Email:** [enquiries@medlinescientific.com](mailto:enquiries@medlinescientific.com)

**Website:** [www.medlinescientific.com](http://www.medlinescientific.com)